



# YZ SERIES

## INFANT OXYGEN HOOD

Stock code [300314]



## Brief Introduction

High transparent acrylic glass material, air & oxygen proportion can be adjusted(YZ-130, YZ-160, YZ-200);  
PC material with beautiful design, resistant to pressure, high strength;  
YZ-300: Two sides operating window can be opened, protect curtain, air valve.

## Product Packaging

Package: The whole unit is packed into one carton  
Dimension: 26cm × 26cm × 28cm (YZ-130&YZ-160)  
32cm × 32cm × 35cm (YZ-200&YZ-300)  
Gross weight: 0.9kg (YZ-130&YZ-160)  
2.0kg (YZ-200&YZ-300)



Ningbo David Medical Device Co.,Ltd

Sales center: #100 Jinghua Road, High-tech Industrial Development Zone, Ningbo, China  
Tel: 0086-574-8780 0008 8780 0007 /Fax: 0086-574-8780 3714  
E-mail: sales@nbdaavid.com /Http: //www.chinadaavid.cn /P.C.: 315040